

PHILIPS

Facial Cleansing
Brush Replacement

SmartClick
accessory

Deep facial cleansing
Brush head replacement



SH560/50



Deep facial cleansing

For smooth and healthy-looking skin

The SmartClick cleansing brush cleans your face thoroughly and leaves your skin with a smooth feeling and a healthy look. So, you feel confident that you always look good!

Smooth and healthy-looking skin


- 32,000 active bristles for thorough cleansing
- Developed for gentle treatment on sensitive skin

Easy to use

- For a hygienic use, replace the brush head every 3 months
- Compatible with SH575 Deep Facial Cleansing Brush

Highlights

32,000 active bristles



The 32,000 active bristles each with a diameter of 50 µm, can reach even the hardest-to-reach areas of your face, cleansing it thoroughly. They remove dirt and dead skin cells, and increase micro-circulation, leaving your skin feeling smooth and looking healthy.

Developed for sensitive skin

The 14mm long bristles have high flexibility, thereby ensuring less friction on sensitive skin, for a gentle cleansing effect. Suitable for everyday use.

Replace head every 3 months

To keep your cleansing brush hygienic, replace the brush head every 3 months by simply clicking on a new brush head replacement (SH560) onto your SmartClick cleansing brush.

Compatible with SH575



Compatible with the SmartClick facial cleansing brushes of Shaver series 7000 and the Deep Facial Cleansing Brush SH575.

Specifications

SmartClick attachment

Fits product type: Shaver series 9000 (S9xxx), Shaver series 7000 (S7xxx), Shaver series 5000 (S5xxx), AquaTouch (S5xxx), Facial cleansing brush SH575

Accessories

Maintenance: Protective cap

Ease of use

Cleaning: Wash with water and mild soap, Dry brush with a towel

Maintenance

Replace brush head: Every 3 months with SH560

