



Philips
In-Ear Headphones

SHE2641



Colour Tunes

Distinctive colour for matching

In-ear headphones with a perfectly balanced, natural sound.

Music to your ears

- Neodymium magnets enhance bass performance and sensitivity
- Twin vents balance the high sounds and bass tones

Made to fit you and your lifestyle

- Ear cushions improve wearing comfort and bass response

Always ready to go

- Durable Flexi-Grip makes a strong flexible connection

PHILIPS

Specifications

Sound

- Acoustic system: open
- Magnet type: Neodymium
- Voice coil: copper
- Diaphragm: Mylar dome
- Frequency response: 12 - 22,000 Hz
- Impedance: 16 ohm
- Maximum power input: 50 mW
- Sensitivity: 103 dB
- Speaker diameter: 13.5 mm
- Type: Dynamic

Connectivity

- Cable Connection: two-parallel, symmetric
- Cable length: 1 m
- Connector: 3.5 mm
- Finishing of connector: chrome-plated
- Type of cable: copper

Inner Carton

- Gross weight: 0.2322 kg
- Inner carton (L x W x H): 22 x 11 x 6.1 cm

- Net weight: 0.063 kg
- Tare weight: 0.1692 kg
- EAN: 87 12581 47862 9
- Number of consumer packages: 6

Outer Carton

- Gross weight: 4.01 kg
- Outer carton (L x W x H): 45.4 x 23.4 x 27.9 cm
- Net weight: 1.008 kg
- Tare weight: 3.002 kg
- EAN: 87 12581 47867 4
- Number of consumer packages: 96

Packaging dimensions

- Packaging type: Blister
- Packaging dimensions (W x H x D): 5 x 17.5 x 2.8 cm
- Gross weight: 0.028 kg
- Net weight: 0.0105 kg
- Tare weight: 0.0175 kg
- EAN: 87 12581 47857 5
- Number of products included: 1

Highlights

Neodymium magnets

Neodymium is the best material for producing a strong magnetic field for greater sensitivity in a voice coil, better bass response and higher overall sound quality.

Twin vents

Specially designed vents that channel the air balance the high tones and bass for a more rounded listening experience.

Comfortable ear cushions

The special shape and luxurious materials used for the ear cushions of these Philips headphones ensure a perfect fit for maximum comfort. They prevent audio leakage and enhance bass performance too. The ear cushions are shaped in such a way that they perfectly align with the area around a person's ear.

Flexi-Grip design

A soft, flexible part protects the cable connection from damage that otherwise would occur due to repeated bending.



Issue date 2013-10-04

Version: 3.0.10

12 NC: 8670 000 51846
EAN: 87 12581 47857 5

© 2013 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com