



Philips
In-Ear Headphones

SHE2642



Colour Tunes

Distinctive colour for matching

In-ear headphones with a perfectly balanced, and natural sound.

Music to your ears

- Neodymium magnets enhance bass performance and sensitivity
- Twin vents balance the high sounds and bass tones

Made to fit you and your lifestyle

- Ear cushions improve wearing comfort and bass response

Always ready to go

- Durable Flexi-Grip makes a strong flexible connection

PHILIPS

Specifications

Sound

- Acoustic system: open
- Magnet type: Neodymium
- Voice coil: copper
- Diaphragm: Mylar dome
- Frequency response: 12 - 22 000 Hz
- Impedance: 16 ohm
- Maximum power input: 50 mW
- Sensitivity: 103 dB
- Speaker diameter: 13.5 mm
- Type: Dynamic

Connectivity

- Cable Connection: two-parallel, symmetric
- Cable length: 1 m
- Connector: 3.5 mm
- Finishing of connector: chrome-plated
- Type of cable: copper

Inner Carton

- Gross weight: 0.193 kg
- Inner carton (L x W x H): 22 x 11 x 6.5 cm

- Nett weight: 0.042 kg
- Tare weight: 0.151 kg
- EAN: 87 12581 47863 6
- Number of consumer packagings: 6

Outer Carton

- Gross weight: 3.596 kg
- Outer carton (L x W x H): 46.5 x 24.5 x 31 cm
- Nett weight: 0.672 kg
- Tare weight: 2.924 kg
- EAN: 87 12581 47868 1
- Number of consumer packagings: 96

Packaging dimensions

- Packaging type: Blister
- Packaging dimensions (W x H x D): 3 x 17.5 x 5.5 cm
- Gross weight: 0.025 kg
- Nett weight: 0.007 kg
- Tare weight: 0.018 kg
- EAN: 87 12581 47858 2
- Number of products included: 1

Highlights

Neodymium magnets

Neodymium is the best material for producing a strong magnetic field for greater sensitivity in a voice coil, better bass response and higher overall sound quality.

Twin vents

Specially designed vents that channel the air balance the high tones and bass for a more rounded listening experience.

Comfort ear cushions

The special shape and luxurious materials used for the ear cushions of these Philips headphones ensure a perfect fit for maximum comfort. They prevent audio leakage and enhance bass performance too. The ear cushions are shaped in such a way that they perfectly align with the area around a person's ear.

Flexi-Grip design

A soft, flexible part protects the cable connection from damage that otherwise would occur due to repeated bending.



Issue date 2014-02-06

Version: 3.0.7

12 NC: 8670 000 51847
EAN: 87 12581 47858 2

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com