



Philips  
In-Ear Headphones

SHE2648



## Colour Tunes

Distinctive colour for matching

In-ear headphones with a perfectly balanced, and natural sound.

### Music to your ears

- Neodymium magnet enhances bass performance and sensitivity
- Twin vents balance the high sounds and bass tones

### Made to fit you and your lifestyle

- Ear cushions improve wearing comfort and bass response

### Always ready to go

- Durable Flexi-Grip makes a strong flexible connection

**PHILIPS**  
sense and simplicity

# Specifications

## Sound

- Acoustic system: open
- Magnet type: Neodymium
- Voice coil: copper
- Diaphragm: Mylar dome
- Frequency response: 12 - 22 000 Hz
- Impedance: 16 ohm
- Maximum power input: 50 mW
- Sensitivity: 103 dB
- Speaker diameter: 13.5 mm
- Type: Dynamic

## Connectivity

- Cable Connection: two-parallel, symmetric
- Cable length: 1 m
- Connector: 3.5 mm
- Finishing of connector: chrome-plated
- Type of cable: copper

## Inner Carton

- Gross weight: 0.2322 kg
- Gross weight: 0.512 lb
- Inner carton (L x W x H): 22 x 11 x 6.1 cm
- Inner carton (L x W x H): 8.7 x 4.3 x 2.4 inch
- Nett weight: 0.139 lb
- Nett weight: 0.063 kg
- Tare weight: 0.1692 kg
- Tare weight: 0.373 lb

- Number of consumer packagings: 6
- GTIN: 2 06 09585 16901 6

## Outer Carton

- Gross weight: 2.407 lb
- Gross weight: 1.092 kg
- Outer carton (L x W x H): 23.3 x 22.8 x 14.5 cm
- Outer carton (L x W x H): 9.2 x 9.0 x 5.7 inch
- Nett weight: 0.556 lb
- Nett weight: 0.252 kg
- Tare weight: 0.84 kg
- Tare weight: 1.852 lb
- Number of consumer packagings: 24
- GTIN: 1 06 09585 16901 9

## Packaging dimensions

- Packaging dimensions (W x H x D): 5 x 17.5 x 2.8 cm
- Packaging dimensions (W x H x D): 2.0 x 6.9 x 1.1 inch
- Gross weight: 0.062 lb
- Gross weight: 0.028 kg
- Nett weight: 0.023 lb
- Nett weight: 0.0105 kg
- Tare weight: 0.0175 kg
- Tare weight: 0.039 lb
- Number of products included: 1
- Packaging type: Blister
- UPC: 6 09585 16901 2

# Highlights

## Neodymium magnet

Neodymium is the best material for producing a strong magnetic field for greater sensitivity in a voice coil, better bass response and higher overall sound quality.

## Twin vents

Specially designed vents that channel the air balance the high tones and bass for a more rounded listening experience.

## Comfort ear cushions

The special shape and luxurious materials used for the ear cushions of these Philips headphones ensure a perfect fit for maximum comfort. They prevent audio leakage and enhance bass performance too. The ear cushions are shaped in such a way that they perfectly align with the area around a person's ear.

## Flexi-Grip design

A soft, flexible part protects the cable connection from damage that otherwise would occur due to repeated bending.



Issue date 2013-01-23

Version: 4.0.5

12 NC: 8670 000 51862  
UPC: 6 09585 16901 2

© 2013 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)