



Balanced sound

for your mobile phone

Twin vents bring greater balance and harmony to your phone's music, while the integrated mic picks up your voice clearly. And with the headsets' soft ear buds you'll enjoy many hours of music and calls in comfort.

Music to your ears

- Neodymium magnet enhances bass performance and sensitivity
- Twin vents balance the high sounds and bass tones

Supreme comfort

- Ear cushions improve wearing comfort and bass response

Designed around you

- Integrated microphone & call button
- 3.5mm plug for Blackberry, HTC, iPhone, LG, Motorola, Palm

Use them for longer

- Durable Flexi-Grip makes a strong flexible connection

Philips
In-Ear Headset



SHE2675OP



PHILIPS

sense and simplicity

Specifications

Sound

- Acoustic system: Open
- Magnet type: Neodymium
- Voice coil: Copper
- Diaphragm: Mylar dome
- Frequency response: 12 - 22 000 Hz
- Impedance: 16 Ohm
- Maximum power input: 50 mW
- Sensitivity: 103 dB
- Speaker diameter: 15 mm
- Type: Dynamic

Connectivity

- Cable Connection: Two-parallel, symmetric
- Cable length: 1.0 m
- Connector: 3.5 mm stereo
- Type of cable: Copper
- Finishing of connector: Chrome plated
- Microphone: Built-in microphone

Packaging dimensions

- Packaging type: Blister
- Number of products included: 1
- Packaging dimensions (W x H x D): 5 x 17 x 3 cm
- Gross weight: 0.0264 kg
- Nett weight: 0.0122 kg
- Tare weight: 0.0142 kg
- Packaging dimensions (W x H x D): 2 x 6.7 x 1.2 inch

- Gross weight: 0.058 lb
- Nett weight: 0.027 lb
- Tare weight: 0.031 lb
- UPC: 6 09585 20665 6

Outer Carton

- Gross weight: 1.2376 kg
- Outer carton (L x W x H): 36.5 x 21.3 x 13.4 cm
- Nett weight: 0.2928 kg
- Number of consumer packagings: 24
- Tare weight: 0.9448 kg
- Gross weight: 2.728 lb
- Outer carton (L x W x H): 14.4 x 8.4 x 5.3 inch
- Nett weight: 0.646 lb
- Tare weight: 2.082 lb
- GTIN: 1 06 09585 20665 3

Inner Carton

- Gross weight: 0.1277 kg
- Inner carton (L x W x H): 17.8 x 9.8 x 5.9 cm
- Nett weight: 0.0366 kg
- Number of consumer packagings: 3
- Tare weight: 0.0911 kg
- Gross weight: 0.282 lb
- Inner carton (L x W x H): 7 x 3.9 x 2.3 inch
- Nett weight: 0.081 lb
- Tare weight: 0.201 lb
- GTIN: 2 06 09585 20665 0

Highlights

Neodymium magnet

Neodymium is the best material for producing a strong magnetic field for greater sensitivity in a voice coil, better bass response and higher overall sound quality.

Twin vents

Specially designed vents that channel the air balance the high tones and bass for a more rounded listening experience.

Comfort ear cushions

The special shape and luxurious materials used for the ear cushions ensure a perfect fit for maximum comfort. They prevent audio leakage and enhance bass performance too. The ear cushions are shaped in such a way that they perfectly align with the area around a person's ear.

Integrated microphone



With the integrated microphone and call button you can use this headset for music as well as calls from your mobile phone. Enjoy handsfree calling, while easily accepting and ending calls from your headset.

3.5mm stereo connector

3.5mm stereo connector for BlackBerry, HTC, iPhone, LG, Motorola and Palm

Flexi-Grip design



A soft, flexible part protects the cable connection from damage that otherwise would occur due to repeated bending.



Issue date 2013-02-09

Version: 3.3.5

12 NC: 8670 000 73005
UPC: 6 09585 20665 6

© 2013 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com