



## Balanced sound

for your mobile phone

Twin vents bring greater balance and harmony to your phone's music, while the integrated mic picks up your voice clearly. And with the headsets' soft ear buds you'll enjoy many hours of music and calls in comfort.

### Music to your ears

- Neodymium magnet enhances bass performance and sensitivity
- Twin vents balance the high sounds and bass tones

### Supreme comfort

- Ear cushions improve wearing comfort and bass response

### Designed around you

- Integrated microphone & call button
- 3.5mm plug for Blackberry, HTC, iPhone, LG, Motorola, Palm

### Use them for longer

- Durable Flexi-Grip makes a strong flexible connection

Philips  
In-Ear Headset



SHE2675YB



# PHILIPS

sense and simplicity

# Specifications

## Sound

- Acoustic system: Open
- Magnet type: Neodymium
- Voice coil: Copper
- Diaphragm: Mylar dome
- Frequency response: 12 - 22 000 Hz
- Impedance: 16 Ohm
- Maximum power input: 50 mW
- Sensitivity: 103 dB
- Speaker diameter: 15 mm
- Type: Dynamic

## Connectivity

- Cable Connection: Two-parallel, symmetric
- Cable length: 1.0 m
- Connector: 3.5 mm stereo
- Type of cable: Copper
- Finishing of connector: Chrome plated
- Microphone: Built-in microphone

## Packaging dimensions

- Number of products included: 1
- Packaging dimensions (W x H x D): 5 x 17 x 3 cm
- Gross weight: 0.0264 kg
- Nett weight: 0.0122 kg
- Tare weight: 0.0142 kg
- Packaging dimensions (W x H x D): 2.0 x 6.7 x 1.2 inch
- Gross weight: 0.058 lb

- Nett weight: 0.027 lb
- Tare weight: 0.031 lb
- UPC: 6 09585 20664 9
- Packaging type: Blister

## Outer Carton

- Gross weight: 1.2376 kg
- Outer carton (L x W x H): 36.5 x 21.3 x 13.4 cm
- Nett weight: 0.2928 kg
- Number of consumer packagings: 24
- Tare weight: 0.9448 kg
- Gross weight: 2.728 lb
- Outer carton (L x W x H): 14.4 x 8.4 x 5.3 inch
- Nett weight: 0.646 lb
- Tare weight: 2.083 lb
- GTIN: 1 06 09585 20664 6

## Inner Carton

- Gross weight: 0.1322 kg
- Inner carton (L x W x H): 17.8 x 9.8 x 5.9 cm
- Nett weight: 0.0366 kg
- Number of consumer packagings: 3
- Tare weight: 0.0956 kg
- Gross weight: 0.291 lb
- Inner carton (L x W x H): 7.0 x 3.9 x 2.3 inch
- Nett weight: 0.081 lb
- Tare weight: 0.211 lb
- GTIN: 2 06 09585 20664 3

# Highlights

## Neodymium magnet

Neodymium is the best material for producing a strong magnetic field for greater sensitivity in a voice coil, better bass response and higher overall sound quality.

## Twin vents

Specially designed vents that channel the air balance the high tones and bass for a more rounded listening experience.

## Comfort ear cushions

The special shape and luxurious materials used for the ear cushions ensure a perfect fit for maximum comfort. They prevent audio leakage and enhance bass performance too. The ear cushions are shaped in such a way that they perfectly align with the area around a person's ear.

## Integrated microphone



With the integrated microphone and call button you can use this headset for music as well as calls from your mobile phone. Enjoy handsfree calling, while easily accepting and ending calls from your headset.

## 3.5mm stereo connector

3.5mm stereo connector for BlackBerry, HTC, iPhone, LG, Motorola and Palm

## Flexi-Grip design



A soft, flexible part protects the cable connection from damage that otherwise would occur due to repeated bending.



Issue date 2012-10-07

Version: 3.2.4

12 NC: 8670 000 73004  
UPC: 6 09585 20664 9

© 2012 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)