



Philips  
In-Ear Headphones

Black



SHE4600



## Enhanced bass

Comfortable fitting with changeable cushioned caps

With 15mm speakers and bass beat vent, these earbud headphones are valued for enhanced bass performance. Airy rubberized soft body provides a comfortable, secure fit inside your ears. 3 sizes of interchangeable earbud for custom fitting.

### Music to your ears

- 15mm speaker driver optimizes wearing comfort
- Bass beat vents allow air movement for better sound
- Neodymium magnets enhance bass performance and sensitivity
- 24k gold-plated plug ensures an ultra reliable connection

### Made to fit you and your lifestyle

- Air cushioned caps for superb cushioning comfort
- Durable Flexi-Grip makes a strong flexible connection

### Always ready to go

- A 1.2m long cable that is ideal for outdoor use

# PHILIPS

# Specifications

## Sound

- Acoustic system: open
- Magnet type: Neodymium
- Diaphragm: Mylar dome
- Frequency response: 18 - 22 000 Hz
- Impedance: 32 ohm
- Maximum power input: 30 mW
- Sensitivity: 107 dB
- Speaker diameter: 15 mm
- Type: Dynamic
- Voice coil: CCAW

## Connectivity

- Cable Connection: two-parallel, symmetric
- Cable length: 1.2 m
- Connector: 3.5 mm
- Finishing of connector: nickel plated
- Type of cable: copper

## Packaging dimensions

- Packaging dimensions (W x H x D): 9,5 x 17,3 x 2,5 cm
- Nett weight: 0,0134 kg
- Gross weight: 0,0592 kg
- Tare weight: 0,0458 kg
- EAN: 69 23410 71659 7
- Number of products included: 1
- Packaging type: Dummy
- Packaging dimensions (W x H x D): 3,7 x 6,8 x 1,0 inch
- Gross weight: 0,131 lb
- Nett weight: 0,030 lb

- Tare weight: 0,101 lb
- Type of shelf placement: Dummy

## Outer Carton

- Outer carton (L x W x H): 35,8 x 21,5 x 18,3 cm
- Nett weight: 0,3216 kg
- Gross weight: 1,99 kg
- Tare weight: 1,6684 kg
- Number of consumer packagings: 24
- GTIN: 1 69 23410 71659 4
- Gross weight: 4,387 lb
- Outer carton (L x W x H): 14,1 x 8,5 x 7,2 inch
- Nett weight: 0,709 lb
- Tare weight: 3,678 lb

## Inner Carton

- Inner carton (L x W x H): 10,5 x 8,8 x 16,9 cm
- Nett weight: 0,0402 kg
- Gross weight: 0,2218 kg
- Tare weight: 0,1816 kg
- Number of consumer packagings: 3
- GTIN: 2 69 23410 71659 1
- Gross weight: 0,489 lb
- Inner carton (L x W x H): 4,1 x 3,5 x 6,7 inch
- Nett weight: 0,089 lb
- Tare weight: 0,400 lb

## Product dimensions

- Product dimensions (W x H x D): 2 x 2,5 x 1,4 cm
- Product dimensions (W x H x D): 0,8 x 1,0 x 0,6 inch
- Weight: 0,0118 kg
- Weight: 0,026 lb

# Highlights

## Air cushioned caps

Air in between the cap and the driver gives you a soft feel like a pillow as it reduces the pressure on your ears.

## 15 mm speaker driver



Small enough for optimum wearing comfort and big enough to deliver crisp, non-distorted sound, the 15 mm speaker driver is an ideal size for listening enjoyment.

## Bass Beat vents

Bass beat vents allow air movement for better sound with a deep rich bass.

## Neodymium magnets

Neodymium is the best material for producing a strong magnetic field for greater sensitivity in a voice coil, better bass response and higher overall sound quality.

## Flexi-Grip design



A soft, flexible part protects the cable connection from damage that otherwise would occur due to repeated bending.

## 24k Gold-plated plug

Rest assured that the precious metal gold finishing on the plug will give you a more reliable connection for better quality audio.

## 1.2m cable

The ideal cable length to give you the freedom to put your audio device where you want.

