

Philips
In-Ear Headphones



SHE9550



Rich bass

Deep, powerful sound

Front acoustics vent provides deep and powerful bass along with an Angled acoustics pipe design which ensures comfortable fit. Included cable winder for hanging onto bags and belts conveniently.

Enjoy quality sound

- Angled Acoustic Pipe channels directional and precise sound
- Exceptional noise isolation for quality sound at low volume
- Acoustics vent control low frequencies for deep clean bass

Ultimate Wearing Comfort

- Ergonomic Angled Acoustic design ensures a natural fit
- 3 sizes of ultra-soft silicon caps ensure an ideal fit
- Soft silicon cushions fit gently in your ears

Always ready to go

- Cable winder simplifies cable storage and length adjustment
- A 1.2m long cable that is ideal for outdoor use
- Durable Flexi-Grip makes a strong flexible connection



PHILIPS

sense and simplicity

Specifications

Sound

- Acoustic system: Semi-open
- Magnet type: Neodymium
- Voice coil: CCAW
- Diaphragm: Mylar dome
- Frequency response: 6 - 23 500 Hz
- Impedance: 16 ohm
- Maximum power input: 50 mW
- Sensitivity: 105 dB
- Speaker diameter: 8.6 mm
- Type: Dynamic

Connectivity

- Cable Connection: symmetrical
- Cable length: 1.2 m
- Connector: 3.5 mm
- Finishing of connector: gold-plated
- Type of cable: OFC

Packaging dimensions

- Type of shelf placement: Dummy
- Packaging dimensions (W x H x D): 9,5 x 17,5 x 3,5 cm
- Nett weight: 0,0395 kg
- Gross weight: 0,083 kg
- Tare weight: 0,0435 kg
- Number of products included: 1
- Packaging type: Dummy
- Packaging dimensions (W x H x D):

- 3,7 x 6,9 x 1,4 inch
- Gross weight: 0,183 lb
- Nett weight: 0,087 lb
- Tare weight: 0,096 lb
- UPC: 6 09585 18711 5

Outer Carton

- Number of consumer packagings: 24
- Gross weight: 6,545 lb
- Gross weight: 2,969 kg
- GTIN: 1 06 09585 18711 2
- Outer carton (L x W x H): 15,3 x 10,1 x 9,8 inch
- Outer carton (L x W x H): 38,8 x 25,6 x 24,8 cm
- Nett weight: 0,9480 kg
- Nett weight: 2,090 lb
- Tare weight: 4,455 lb
- Tare weight: 2,0210 kg

Inner Carton

- Number of consumer packagings: 3
- Gross weight: 0,634 lb
- Gross weight: 0,2876 kg
- GTIN: 2 06 09585 18711 9
- Inner carton (L x W x H): 7,2 x 4,6 x 4,1 inch
- Inner carton (L x W x H): 18,3 x 11,7 x 10,4 cm
- Nett weight: 0,1185 kg
- Nett weight: 0,261 lb
- Tare weight: 0,373 lb
- Tare weight: 0,1691 kg



Issue date 2013-01-22

Version: 1.0.3

12 NC: 8670 000 61019
UPC: 6 09585 18711 5

© 2013 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com