

celia  
miller

Philips O'Neill  
Headband headphones

- O'Neill
- On-ear
- Blue



SHO4300BL

## Incredibly Addictive Bass

Deep Bass. Flat cans.

Enjoy deeper, richer bass notes and clearer vocals with the Snug's 40mm drivers. An ultra-low profile design with incredibly comfy ear cushions that pump out bass makes these flat folding cans just perfect for the longest journeys.

### Supreme comfort

- Comfortable headband easily adjusts to every head size
- Headband's soft inner cushion evens out wearing pressure
- Soft leather ear cushions prolong wearing comfort

### Mind-blowing sound

- 40mm neodymium drivers deliver deeper, more dynamic bass
- Fill your head with rock via on-ear isolating cushions

### Prepare for impact

- Flat, foldable design simplifies storage and portability
- C-hinge design for more comfort and flexibility
- L-shape connector enhances cable durability and lifespan
- Increased strain relief length enables more severe bending

PHILIPS

O'NEILL

Headband headphones  
O'Neill On-ear, Blue

SHO4300BL/00

# Highlights

## 40mm neodymium drivers



Rock out to high quality tunes on the go with dynamic 40mm drivers that deliver crystal clear sound coupled with thumping bass.

## Comfy and adjustable headband



Comfortable headband easily adjusts to every head size

## Flat foldable earcans



Flat, foldable design simplifies storage and portability

## Flexible C-hinge



C-hinge design for more comfort and flexibility

## L-shape connector



Enjoy music for longer with a reinforced L-shaped cable connector that makes these headphones more durable and longer lasting.

## On-ear isolating cushions



Fill your head with rock via on-ear isolating cushions. Not only will you enjoy long lasting comfort, these cushions help you hear and enjoy every sonic detail of your favorite tunes.

## Reinforced strain relief

Increased strain relief length enables more severe bending

## Soft inner headband cushion



Headband's soft inner cushion evens out wearing pressure

## Soft leather ear cushions



Ultra-soft leather ear cushions hug the contours of your ears comfortably and securely, creating a perfect seal that is so comfortable you might forget you have headphones on.

# Specifications

## Sound

- Acoustic system: closed
- Magnet type: Neodymium
- Voice coil: Copper
- Diaphragm: PET
- Impedance: 32 ohm
- Maximum power input: 50 mW
- Sensitivity: 105 dB
- Speaker diameter: 40 mm
- Type: Dynamic
- Frequency response: 18 - 22 000 Hz

## Connectivity

- Cable Connection: single-sided cable
- Cable length: 1.2 m
- Connector: 3.5 mm
- Finishing of connector: gold-plated
- Type of cable: OFC

## Outer Carton

- Gross weight: 2.511 lb
- Gross weight: 1.139 kg
- GTIN: 1 69 23410 72090 4
- Outer carton (L x W x H): 26 x 19 x 22.5 cm
- Outer carton (L x W x H): 10.2 x 7.5 x 8.9 inch
- Nett weight: 0.860 lb
- Nett weight: 0.39 kg

- Number of consumer packagings: 3
- Tare weight: 1.651 lb
- Tare weight: 0.749 kg

## Packaging dimensions

- Packaging dimensions (W x H x D):  
17 x 22.2 x 8 cm
- Packaging dimensions (W x H x D):  
6.7 x 8.7 x 3.1 inch
- EAN: 69 23410 72090 7
- Gross weight: 0.283 kg
- Gross weight: 0.624 lb
- Nett weight: 0.287 lb
- Nett weight: 0.13 kg
- Number of products included: 1
- Packaging type: Carton
- Tare weight: 0.153 kg
- Tare weight: 0.337 lb
- Type of shelf placement: Both

## Product dimensions

- Product dimensions (W x H x D):  
5.9 x 7.4 x 1.5 inch
- Product dimensions (W x H x D):  
15.1 x 18.7 x 3.8 cm
- Weight: 0.122 kg
- Weight: 0.269 lb