



Philips O'Neill
Headband headphones

- O'Neill
- On-ear
- Pink

SHO4300PK



Incredibly Addictive Bass

Deep Bass. Flat cans.

Enjoy deeper, richer bass notes and clearer vocals with the Snug's 40mm drivers. An ultra-low profile design with incredibly comfy ear cushions that pump out bass makes these flat folding cans just perfect for the longest journeys.

Supreme comfort

- Comfortable headband easily adjusts to every head size
- Headband's soft inner cushion evens out wearing pressure
- Soft leather ear cushions prolong wearing comfort

Mind-blowing sound

- 40mm neodymium drivers deliver deeper, more dynamic bass
- Fill your head with rock via on-ear isolating cushions

Prepare for impact

- Flat, foldable design simplifies storage and portability
- C-hinge design for more comfort and flexibility
- L-shape connector enhances cable durability and lifespan
- Increased strain relief length enables more severe bending

PHILIPS

O'NEILL

Headband headphones
O'Neill On-ear, Pink

SHO4300PK/00

Highlights

40mm neodymium drivers



Rock out to high quality tunes on the go with dynamic 40mm drivers that deliver crystal clear sound coupled with thumping bass.

Comfy and adjustable headband



Comfortable headband easily adjusts to every head size

Flat foldable earcans



Flat, foldable design simplifies storage and portability

Flexible C-hinge



C-hinge design for more comfort and flexibility

L-shape connector



Enjoy music for longer with a reinforced L-shaped cable connector that makes these headphones more durable and longer lasting.

On-ear isolating cushions



Fill your head with rock via on-ear isolating cushions. Not only will you enjoy long lasting comfort, these cushions help you hear and enjoy every sonic detail of your favorite tunes.

Reinforced strain relief

Increased strain relief length enables more severe bending

Soft inner headband cushion



Headband's soft inner cushion evens out wearing pressure

Soft leather ear cushions



Ultra-soft leather ear cushions hug the contours of your ears comfortably and securely, creating a perfect seal that is so comfortable you might forget you have headphones on.

Specifications

Sound

- Acoustic system: closed
- Magnet type: Neodymium
- Voice coil: Copper
- Diaphragm: PET
- Impedance: 32 ohm
- Maximum power input: 50 mW
- Sensitivity: 105 dB
- Speaker diameter: 40 mm
- Type: Dynamic
- Frequency response: 18 - 22 000 Hz

Connectivity

- Cable Connection: single-sided cable
- Cable length: 1.2 m
- Connector: 3.5 mm
- Finishing of connector: gold-plated
- Type of cable: OFC

Outer Carton

- Gross weight: 2.511 lb
- Gross weight: 1.139 kg
- GTIN: 1 69 23410 72089 8
- Outer carton (L x W x H): 26 x 19 x 22.5 cm
- Outer carton (L x W x H): 10.2 x 7.5 x 8.9 inch
- Nett weight: 0.860 lb
- Nett weight: 0.39 kg

- Number of consumer packagings: 3
- Tare weight: 1.651 lb
- Tare weight: 0.749 kg

Packaging dimensions

- Packaging dimensions (W x H x D):
17 x 22.2 x 8 cm
- Packaging dimensions (W x H x D):
6.7 x 8.7 x 3.1 inch
- EAN: 69 23410 72089 1
- Gross weight: 0.283 kg
- Gross weight: 0.624 lb
- Nett weight: 0.287 lb
- Nett weight: 0.13 kg
- Number of products included: 1
- Packaging type: Carton
- Tare weight: 0.153 kg
- Tare weight: 0.337 lb
- Type of shelf placement: Both

Product dimensions

- Product dimensions (W x H x D):
5.9 x 7.4 x 1.5 inch
- Product dimensions (W x H x D):
15.1 x 18.7 x 3.8 cm
- Weight: 0.122 kg
- Weight: 0.269 lb