




Philips O'Neill
THE SNUG headband
headphones

The Snug

Flawless fit and sound

The flat-folding Snug is perfect for surf trips: just fold 'em & throw 'em into your bag. Ultrasoft cushioning & next level engineering team up to produce an acoustic experience of the highest order.



SHO8803

Powerful sound

- Dynamic 40mm drivers deliver deep, dynamic bass

Enhanced durability

- Cable stress relief to prevent breakage
- Tough tangle-free cable and reinforced connectors

Supreme comfort


- Sound isolating super soft ear cushions
- Adjustable fit
- Folds flat for easy storing and carrying

PHILIPS

O'NEILL


Highlights

Smart reinforcement




Thanks to the tough tangle-free fabric cable and reinforced connectors these headphones are built to stand up to the wear and tear of your active lifestyle.

Folds Flat for portability



Folds flat for easy storing and carrying

Cable stress relief



Inspired by the latest in ski technology the cable quickly disconnects from the headband whenever

the cable is caught. Simply reconnect the cable to the cable stress relief of the headband and you can continue enjoying your music.

Powerful sound

With dynamic 40mm drivers these headphones deliver a crystal clear sound experience coupled with a bass thumping performance.

Comfortable for long wear



The deluxe super soft ear cushions enable you to wear your headphones for as long as you like without discomfort.

Adjustable fit



The headband of these Philips headphones can easily be adjusted to snugly fit any head any time.

Specifications

Sound

- Acoustic system: Closed
- Magnet type: Neodymium
- Voice coil: CCAW
- Diaphragm: Mylar dome
- Frequency response: 15 - 24 000 Hz
- Impedance: 32 Ohm
- Maximum power input: 100 mW
- Sensitivity: 107 dB
- Speaker diameter: 40 mm
- Type: Dynamic

Connectivity

- Cable Connection: Single-sided cable
- Cable length: 1.2 m
- Connector: 3.5 mm stereo
- Finishing of connector: 24k Gold plated
- Type of cable: OFC

Outer Carton

- Gross weight: 1.175 kg
- Gross weight: 2.590 lb
- GTIN: 1 06 09585 19453 0
- Outer carton (L x W x H): 8.3 x 6.6 x 9.4 inch
- Outer carton (L x W x H): 21 x 16.8 x 23.8 cm
- Nett weight: 0.444 kg
- Nett weight: 0.979 lb
- Number of consumer packagings: 3
- Tare weight: 0.731 kg
- Tare weight: 1.612 lb

Packaging dimensions

- Packaging dimensions (W x H x D): 19.5 x 21.5 x 5 cm
- Packaging dimensions (W x H x D): 7.7 x 8.5 x 2.0 inch
- Gross weight: 0.32 kg
- Gross weight: 0.705 lb
- Nett weight: 0.326 lb
- Nett weight: 0.148 kg
- Number of products included: 1
- Packaging type: Blister
- Tare weight: 0.379 lb
- Tare weight: 0.172 kg
- UPC: 6 09585 19453 3

