




Philips ActionFit
Sports earhook
headphones

ActionFit
Orange & Grey

SHQ3005



Tuned for sports

Adjustable fit

Constructed with water proof materials coupled with an anti bacterial agent that helps kill germs, these headphones have been designed for the fitness enthusiast

Innovation for Athletes

- Customizable ear hooks
- Ultra-soft silicon caps sit comfortably in your ears
- A 1.2m long cable that is ideal for outdoor use

Designed for active use

- Washable
- Sweat proof
- Cable clip

Energize you

- Bass thumping stereo sound - fuel for your workout
- Sound Isolation headphones to minimize environmental noise

PHILIPS

Highlights

Washable



Constructed with water proof materials allowing you to wash your headphones under running tap water

Sweat proof



Designed with a water resistant seal to protect your headphones from sweat related damage

Customizable ear hooks



Customizable ear hooks for enhanced fit - designed for sports

Bass thumping stereo sound

Acoustics finely tuned to enhance the sporting experience - fuel for your workout

3-size ultra-soft silicon caps



The 3 sizes of ultra-soft silicon caps of these Philips earphones sit comfortably in your ears without adding pressure to them whilst blocking out unwanted noise.

Cable clip



Clip helps fix cable securely

Sound Isolation



With sound isolation, you will hear and enjoy greater music details at a lower volume than is possible with normal headphones. Based on thorough understanding of the human ear canal anatomy, the sound isolation earphone's angled acoustic design creates a perfect seal between the ambient sound source and the ear canal, thus eliminating unwanted background noise while delivering sound directly to your ears. Even at a moderate volume setting, you can still discern music details without noise distractions from the environment. Earphone sleeves are available in three sizes to provide the best fit.

1.2m cable

The ideal cable length to give you the freedom to put your audio device where you want.

Specifications

Sound

- Diaphragm: Mylar dome
- Frequency response: 15 - 22 000 Hz
- Impedance: 16 ohm
- Magnet type: Neodymium
- Maximum power input: 20 mW
- Sensitivity: 102 dB
- Speaker diameter: 9 mm
- Type: Dynamic
- Voice coil: CCAW

Connectivity

- Cable Connection: two-parallel, symmetric
- Cable length: 1.2 m
- Connector: 3.5 mm
- Finishing of connector: gold-plated
- Type of cable: OFC

Accessories

- Included: transport case

Packaging dimensions

- Packaging dimensions (W x H x D):
9.5 x 17.5 x 3.5 cm
- Nett weight: 0.024 kg
- Gross weight: 0.075 kg

- Tare weight: 0.051 kg
- EAN: 87 12581 63853 5
- Number of products included: 1
- Packaging type: Carton
- Type of shelf placement: Both

Outer Carton

- Outer carton (L x W x H): 38.8 x 25.6 x 24.8 cm
- Nett weight: 0.576 kg
- Gross weight: 2.587 kg
- Tare weight: 2.011 kg
- Number of consumer packagings: 24
- GTIN: 1 87 12581 63853 2

Inner Carton

- Inner carton (L x W x H): 18.3 x 11.7 x 10.4 cm
- Nett weight: 0.072 kg
- Gross weight: 0.27 kg
- Tare weight: 0.198 kg
- Number of consumer packagings: 3
- GTIN: 2 87 12581 63853 9

Product dimensions

- Product dimensions (W x H x D):
3.95 x 6.35 x 2.4 cm
- Weight: 0.024 kg



Issue date 2014-03-17

Version: 1.1.4

12 NC: 8670 000 88116
EAN: 87 12581 63853 5

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com