

Philips
Earhook Headphones



SHS4841



Extra comfortable

Adjustable earhook fits all ears

This inner earhook sleeves of this headphone can be adjusted to sit on 5 different positions for a more personalized fit. Air in-between the caps and speakers allows for a soft feel as the caps become softer like a cushion.

Music to your ears

- 24k gold-plated plug ensures an ultra reliable connection
- Bass beat vents allow air movement for better sound

Made to fit you and your lifestyle

- Contoured earhook increase wearing comfort and stability
- 13.5mm speaker driver optimizes wearing comfort
- Air cushioned caps for superb cushioning comfort

Always ready to go

- Stays comfortably in your ear when you are on the move
- Durable Flexi-Grip makes a strong flexible connection

PHILIPS

Specifications

Product dimensions

- Product dimensions (W x H x D): 4,2 x 6,8 x 1,7 cm
- Product dimensions (W x H x D): 1,7 x 2,7 x 0,7 inch
- Weight: 0,0151 kg
- Weight: 0,033 lb

Sound

- Acoustic system: open
- Magnet type: Neodymium
- Diaphragm: Mylar dome
- Voice coil: copper
- Frequency response: 12 - 23 500 Hz
- Sensitivity: 105 dB
- Impedance: 16 ohm
- Maximum power input: 50 mW
- Speaker diameter: 13.5 mm
- Type: Dynamic

Connectivity

- Cable Connection: two-parallel, symmetric
- Cable length: 1.2 m
- Connector: 3.5 mm
- Finishing of connector: chrome-plated
- Type of cable: OFC

Inner Carton

- Gross weight: 0,2886 kg
- Gross weight: 0,637 lb
- GTIN: 2 06 09585 18004 2
- Inner carton (L x W x H): 8,1 x 7,0 x 2,5 inch
- Inner carton (L x W x H): 20,6 x 17,9 x 6,3 cm

- Nett weight: 0,0906 kg
- Nett weight: 0,200 lb
- Number of consumer packagings: 6
- Tare weight: 0,437 lb
- Tare weight: 0,1980 kg

Outer Carton

- Gross weight: 1,2944 kg
- Gross weight: 2,854 lb
- GTIN: 1 06 09585 18004 5
- Outer carton (L x W x H): 14,5 x 8,4 x 5,5 inch
- Outer carton (L x W x H): 36,8 x 21,4 x 14 cm
- Nett weight: 0,3624 kg
- Nett weight: 0,799 lb
- Number of consumer packagings: 24
- Tare weight: 2,055 lb
- Tare weight: 0,9320 kg

Packaging dimensions

- Type of shelf placement: Dummy
- Packaging dimensions (W x H x D): 5,3 x 17,3 x 3,8 cm
- Gross weight: 0,0356 kg
- Nett weight: 0,0151 kg
- Tare weight: 0,0205 kg
- Packaging dimensions (W x H x D): 2,1 x 6,8 x 1,5 inch
- Gross weight: 0,078 lb
- Nett weight: 0,033 lb
- Number of products included: 1
- Packaging type: Dummy
- Tare weight: 0,045 lb
- UPC: 6 09585 18004 8

Highlights

Contoured earhook

The earhook stays put behind your ear even when you're on the move, keeping the earphone comfortably in place and sounding great. It's the part contacting the back of your ears and keeps the headphones on your ears.

Designed for active use

Whether you are out running, cycling or taking part in any active pastime, this headphone will stay comfortably in your ear.



Issue date 2014-08-21

Version: 2.0.10

12 NC: 8670 000 57991
UPC: 6 09585 18004 8

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com