

Philips
Earhook Headphones



SHS4842



Extra comfortable

Adjustable earhook fits all ears

This inner earhook sleeves of this headphone can be adjusted to sit on 5 different positions for a more personalized fit. Air in-between the caps and speakers allows for a soft feel as the caps become softer like a cushion.

Music to your ears

- 24k gold-plated plug ensures an ultra reliable connection
- Bass beat vents allow air movement for better sound

Made to fit you and your lifestyle

- Contoured earhook increase wearing comfort and stability
- 13.5mm speaker driver optimizes wearing comfort
- Air cushioned caps for superb cushioning comfort

Always ready to go

- Stays comfortably in your ear when you are on the move
- Durable Flexi-Grip makes a strong flexible connection

PHILIPS
sense and simplicity

Specifications

Product dimensions

- Product dimensions (W x H x D): 3,6 x 6 x 1,8 cm
- Product dimensions (W x H x D): 1,4 x 2,4 x 0,7 inch
- Weight: 0,016 kg
- Weight: 0,035 lb

Sound

- Acoustic system: open
- Magnet type: Neodymium
- Diaphragm: Mylar dome
- Voice coil: copper
- Frequency response: 12 - 23 500 Hz
- Sensitivity: 105 dB
- Impedance: 16 ohm
- Speaker diameter: 13.5 mm
- Type: Dynamic
- Maximum power input: 15 mW

Connectivity

- Cable Connection: two-parallel, symmetric
- Cable length: 1.2 m
- Connector: 3.5 mm
- Finishing of connector: chrome-plated
- Type of cable: OFC

Inner Carton

- GTIN: 2 06 09585 18735 5
- Inner carton (L x W x H): 7,8 x 7,3 x 2,6 inch
- Inner carton (L x W x H): 19,8 x 18,5 x 6,7 cm
- Number of consumer packagings: 6
- Gross weight: 0,697 lb

- Gross weight: 0,316 kg
- Nett weight: 0,0966 kg
- Nett weight: 0,213 lb
- Tare weight: 0,2194 kg
- Tare weight: 0,484 lb

Outer Carton

- Gross weight: 3,448 lb
- Gross weight: 1,564 kg
- GTIN: 1 06 09585 18735 8
- Outer carton (L x W x H): 15,4 x 8,7 x 6,9 inch
- Outer carton (L x W x H): 39,2 x 22 x 17,4 cm
- Nett weight: 0,852 lb
- Nett weight: 0,3864 kg
- Number of consumer packagings: 24
- Tare weight: 2,596 lb
- Tare weight: 1,1776 kg

Packaging dimensions

- Type of shelf placement: Dummy
- Packaging dimensions (W x H x D): 5,2 x 17 x 3 cm
- Gross weight: 0,036 kg
- Nett weight: 0,0161 kg
- Tare weight: 0,0199 kg
- Packaging dimensions (W x H x D): 2,0 x 6,7 x 1,2 inch
- Gross weight: 0,079 lb
- Nett weight: 0,035 lb
- Number of products included: 1
- Packaging type: Dummy
- Tare weight: 0,044 lb
- UPC: 6 09585 18735 1

Highlights

Contoured earhook

The earhook stays put behind your ear even when you're on the move, keeping the earphone comfortably in place and sounding great. It's the part contacting the back of your ears and keep the headphones on your ears.

Designed for active use

Whether you are out running, cycling or taking part in any active pastime, this headphone will stay comfortably in your ear.

13.5mm speaker driver

Small enough for optimum wearing comfort and big enough to deliver undistorted sound, an ideal size for listening enjoyment.

24k Gold-plated plug

Rest assured that the precious metal gold finishing on the plug will give you a more reliable connection for better quality audio.


Air cushioned caps

Air in between the cap and the driver gives you a soft feel like a pillow as it reduces the pressure on your ears.

Bass beat vents

Bass beat vents allow air movement for better sound with a deep rich bass.

Flexi-grip design



The soft flexi-grip material that the body of the headphone is sheathed in enhances its performance and extends its durability. It protects the cable connection from damage that otherwise could occur due to repeated bending of the headphone body from frequent contact with the user's fingertips when being put on and removed.



Issue date 2013-01-22

Version: 2.0.7

12 NC: 8670 000 61126
UPC: 6 09585 18735 1

© 2013 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com