



Philips
Leather wallet

For iPod nano



SJM3202



Protect your nano in style

Stylish leather fold-over design protects your nano and conveniently stores cash and cards.

Protect your iPod

- Protective leather design guards against scratches

PHILIPS

Leather wallet
For iPod nano

SJM3202/10

Specifications

Packaging

- EAN/UPC/GTIN: 87 10895 98647 2
- Length: 24 cm
- Width: 11.2 cm
- Height: 2.3 cm
- Gross Weight: .04 kg
- Tare Weight: .02 kg
- Net Weight: .02 kg

Inner Carton

- EAN/UPC/GTIN: 1 87 10895 98647 9
- Quantity: 4
- Length (cm): 25.5 cm
- Width (cm): 12 cm
- Height (cm): 12 cm
- Gross Weight: .26 kg
- Tare Weight: .08 kg
- Net Weight: .18 kg

Outer Carton

- EAN/UPC/GTIN: 2 87 10895 98647 6
- Quantity: 24
- Length (cm): 38 cm
- Width (cm): 26.5 cm

- Height (cm): 27.5 cm
- Gross Weight: 2.08 kg
- Tare Weight: .48 kg
- Net Weight: 1.6 kg

Packaging dimensions

- Packaging dimensions (W x H x D):
11.5 x 25.5 x 3.5 cm
- Gross weight: .1141 kg
- Tare weight: .0941 kg
- Nett weight: .02 kg

Outer Carton

- Outer carton (L x W x H): 40 x 29.5 x 28.5 cm
- Gross weight: 2.74 kg
- Tare weight: 2.26 kg
- Nett weight: .48 kg

Inner Carton

- Inner carton (L x W x H): 13.5 x 12.5 x 12 cm
- Gross weight: .4567 kg
- Tare weight: .3767 kg
- Nett weight: .08 kg

Highlights

Protective leather design

This lightweight leather case protects against scratches and normal wear and tear.



Issue date 2014-02-06

Version: 2.1.8

12 NC: 9082 100 09263
EAN: 87 10895 98647 2

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com