



Philips
Screen cleaner

Notebook

SPC2500

Keep your notebook screen clean

Clean and maintain your LCD flat screen and laptop monitors. Uses an anti-static cleaning gel and soft cleaning cloth that is safe and effective.

Safely clean your screen

- Streak-free cleaning gel for spotless screens
- Lint-free microfiber cloth cleans safely and thoroughly

PHILIPS

Specifications

Packaging

- EAN/UPC/GTIN: 87 12581 30613 7
- Quantity: 1
- Gross Weight: 0.110 kg
- Tare Weight: 0.020 kg
- Net Weight: 0.090 kg
- Length: 2.8 cm
- Width: 13.1 cm
- Height: 19.0 cm

Outer Carton

- EAN/UPC/GTIN: 87 12581 30697 7
- Quantity: 6
- Gross weight: 0.820 kg
- Tare weight: 0.100 kg
- Net weight: 0.720 kg
- Length (cm): 18.0 cm
- Width (cm): 15.0 cm
- Height (cm): 22.0 cm

Highlights

Streak-free cleaning gel

Enjoy a spotless screen with this streak-free cleaning gel, specially formulated not to leave any kind of residue on your screen.

Lint-free microfiber cloth

Safely and thoroughly clean your screen with this lint-free microfiber cloth, designed to eliminate dust, dirt and fingerprints without leaving anything behind.



Issue date 2014-02-07

Version: 3.1.8

12 NC: 8670 000 25911
EAN: 87 12581 30613 7

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com