

Philips
Wired notebook mouse

USB
800 DPI

SPM1702VB



Wired notebook mouse

Retractable cable

The retractable cable of this notebook mouse offers you convenience wherever you are. The accuracy is equally impressive, with 800DPI and a high definition precision optical sensor guaranteeing you smooth and accurate tracking.

Powerful performance

- 800 DPI accurate reactions
- High-performance LED optic guarantees precision control

Excellent ease-of-use

- This mouse with retractable cable is ideal for travellers

PHILIPS

Wired notebook mouse
USB 800 DPI

SPM1702VB/97

Specifications

System Requirements

- Operating System: Microsoft Windows 2000, Microsoft Windows ME, Microsoft Windows XP, Microsoft Windows Vista

Connectivity

- PC Link: USB interface

Accessories

- Included accessories: IFU / User Manual

Product dimensions

- Product dimensions (W x H x D): 4.9 x 9 x 3.2 cm
- Weight: .055 kg

Packaging dimensions

- Packaging dimensions (W x H x D): 16.9 x 23.5 x 4.2 cm
- Nett weight: .055 kg
- Gross weight: .15 kg
- Tare weight: .095 kg

Outer Carton

- Outer carton (L x W x H): 51 x 47 x 20 cm
- Nett weight: 1.1 kg
- Gross weight: 4.05 kg
- Tare weight: 2.95 kg

Highlights

800 DPI accurate reactions

With 800DPI, this mouse provides accurate reactions. The cursor will move exactly where you want it. It's all thanks to the high-definition optical sensor that registers even the smallest movements. Giving you complete cursor control.

Precise optical sensor

Philips optical technology ensures extra accuracy and responsiveness. You can work with it on most surfaces, including wood, plastic or paper.

Retractable cable mouse

This retractable cable saves space making it ideal for travellers



Issue date 2014-02-05

Version: 1.0.8

12 NC: 8670 000 34239
EAN: 87 12581 36632 2

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com