



Philips Cinema 21:9
TV stand

CINEMA^{21:9}

ST569954

Comfortable viewing of your Cinema 21:9 TV

This high quality stand in matching glass and metal finish allows optimal placement in your interior

Freedom of placement and use

- Convenient 30° swivel
- 2-in-1 stand for easy use on table and as swivel wall mount

Stylish and elegant design

- Elegant yet strong metal support
- High-quality glass top plate
- Matching design with Philips Cinema 21:9

PHILIPS
sense and simplicity

Specifications

Related Products

- Compatible with: 56PFL9954

Design

- Colour: Black

Dimensions

- Stand dimensions (W x H x D):

- 565 x 200 x 325 mm
- Box dimensions (W x H x D):
660 x 280 x 390 mm
- Weight incl. Packaging: 8.5 kg

Packaging Data

- Delivered as: Fully assembled



Issue date 2011-02-02

Version: 1.2.5

12 NC: 8670 000 51303
EAN: 87 12581 47541 3

© 2011 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com