



Philips
Storage cases

Clear

STO3500

CD & DVD protective sleeves

Safeguard your media. Double-sided design for twice the storage. Quantity 50.

Protect your CDs and DVDs

- Sturdy construction

Store more discs

- Double storage design

PHILIPS

Storage cases
Clear

STO3500/97

Specifications

Inner Carton

- Quantity: 6
- Gross weight: 1.68 kg
- Tare weight: 0.126 kg
- Net weight: 1.55 kg
- EAN/UPC/GTIN: 87 12581 35254 7
- Height (cm): 17.5 cm
- Length (cm): 24.5 cm
- Width (cm): 16.5 cm

Outer Carton

- Quantity: 72
- Gross weight: 21.4 kg
- Tare weight: 2.75 kg
- Net weight: 18.64 kg

- EAN/UPC/GTIN: 87 12581 35255 4
- Height (cm): 54.5 cm
- Length (cm): 51.0 cm
- Width (cm): 34.5 cm

Packaging

- EAN/UPC/GTIN: 87 12581 35253 0
- Gross Weight: 0.260 kg
- Net Weight: 0.259 kg
- Tare Weight: 0.001 kg
- Height: 17.0 cm
- Width: 15.0 cm
- Length: 4.0 cm
- Quantity: 1

Highlights

Sturdy construction

Strong plastic transparent cases store and protect your discs from dust and dirt.

Double storage design

Since a second side allows you to store a second disc in the sleeve, you can store more discs in the same amount of space.



Issue date 2014-02-07

Version: 2.0.3

12 NC: 8670 000 32522
EAN: 87 12581 35253 0

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com