


Philips  
Storage cases



STO3520



## CD & DVD protective sleeves

Safeguard your media on the go. Different color sleeves make it easy to identify discs.  
Quantity 50.

### **Protect your CDs and DVDs**

- Protective paper sleeve

### **Store more discs**

- Space saving design

**PHILIPS**  
sense and simplicity

# Specifications

## Packaging Data

- Quantity: 1
- EAN/UPC/GTIN: 6 09585 13311 2
- Length: 160 mm
- Width: 130 mm
- Height: 40 mm
- Gross weight: 0.16 kg
- Tare weight: .002 kg
- Net weight: 1.58 kg

## Inner Carton

- EAN/UPC/GTIN: 1 06 09585 13311 9
- Quantity: 6
- Length: 215 mm
- Width: 145 mm

- Height: 165 mm
- Gross weight: 1.14 kg
- Tare weight: 0.100 kg
- Net weight: 1.040 kg

## Outer Carton

- EAN/UPC/GTIN: 2 06 09585 13311 6
- Quantity: 72
- Length: 450 mm
- Width: 305 mm
- Height: 525 mm
- Gross weight: 14.6 kg
- Tare weight: 0.200 kg
- Net weight: 144 kg

# Highlights

## Protective sleeve

Paper sleeve shields the disc from dust and contaminants

## Space saving design

Since paper sleeves are thinner than jewel cases, you can store more discs in the same amount of space.



Issue date 2009-11-06

Version: 1.1.8

12 NC: 8670 000 30023  
UPC: 6 09585 13311 2

© 2009 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips  
Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)