

Philips
Storage cases



STO3520

CD & DVD protective sleeves

Safeguard your media on the go. Different color sleeves make it easy to identify discs.
Quantity 50.

Protect your CDs and DVDs

- Protective paper sleeve

Store more discs

- Space saving design

PHILIPS

Specifications

Inner Carton

- Quantity: 6
- Gross weight: 1.14 kg
- Tare weight: 0.100 kg
- EAN/UPC/GTIN: 87 12581 35257 8
- Height (cm): 16.5 cm
- Length (cm): 21.5 cm
- Net weight: 1.038 kg
- Width (cm): 14.0 cm

Outer Carton

- Quantity: 72
- Gross weight: 14.6 kg
- Tare weight: 2.2 kg
- Net weight: 12.5 kg

- EAN/UPC/GTIN: 87 12581 35258 5
- Height (cm): 52.5 cm
- Length (cm): 45.0 cm
- Width (cm): 30.5 cm

Packaging

- EAN/UPC/GTIN: 87 12581 35256 1
- Gross Weight: 0.18 kg
- Height: 16.0 cm
- Length: 13.0 cm
- Net Weight: 0.173 kg
- Quantity: 1
- Tare Weight: 0.007 kg
- Width: 4.0 cm

Highlights

Protective sleeve

Paper sleeve shields the disc from dust and contaminants

Space saving design

Since paper sleeves are thinner than jewel cases, you can store more discs in the same amount of space.



Issue date 2014-02-07

Version: 2.0.4

12 NC: 8670 000 32523
EAN: 87 12581 35256 1

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com