

Philips  
Storage cases



STO3520

## CD & DVD protective sleeves

Safeguard your media on the go. Different color sleeves make it easy to identify discs.  
Quantity 50.

### **Protect your CDs and DVDs**

- Protective paper sleeve

### **Store more discs**

- Space saving design

**PHILIPS**  
sense and simplicity

# Specifications

## Inner Carton

- Quantity: 6
- Gross weight: 1.14 kg
- Tare weight: 0.100 kg
- EAN/UPC/GTIN: 87 12581 35257 8
- Height (cm): 16.5 cm
- Length (cm): 21.5 cm
- Net weight: 1.038 kg
- Width (cm): 14.0 cm

## Outer Carton

- Quantity: 72
- Gross weight: 14.6 kg
- Tare weight: 2.2 kg
- Net weight: 12.5 kg

- EAN/UPC/GTIN: 87 12581 35258 5
- Height (cm): 52.5 cm
- Length (cm): 45.0 cm
- Width (cm): 30.5 cm

## Packaging

- EAN/UPC/GTIN: 87 12581 35256 1
- Gross Weight: 0.18 kg
- Height: 16.0 cm
- Length: 13.0 cm
- Net Weight: 0.173 kg
- Quantity: 1
- Tare Weight: 0.007 kg
- Width: 4.0 cm

# Highlights

## Protective sleeve

Paper sleeve shields the disc from dust and contaminants

## Space saving design

Since paper sleeves are thinner than jewel cases, you can store more discs in the same amount of space.



Issue date 2009-11-06

Version: 2.0.4

12 NC: 8670 000 32523  
EAN: 87 12581 35256 1

© 2009 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips  
Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)