



Bass Sound

for enhanced sound

With enhanced bass performance, these in-ear headphones bring quality music comfortably to your ears.

Enjoy quality sound

- Neodymium speaker drivers deliver pure balanced sound
- 14.2mm speaker drivers for rich bass and clear sound

Ultimate wearing comfort

- This ergonomic design will fit perfectly into your ear canal
- Designed for ear geometry comfort fit

Philips
Headphones

Earbud



TAUE100BK



PHILIPS

Specifications

Design

- Color: Black

Sound

- Impedance: 32 ohm
- Maximum power input: 5 mW
- Speaker diameter: 14.2 mm
- Sensitivity: 106 dB
- Frequency response: 20 - 20 000 Hz
- Acoustic system: open

Connectivity

- Cable length: 1 m

Inner Carton

- Gross weight: 0.402 kg
- GTIN: 2 48 95229 10055 2
- Inner carton (L x W x H): 31 x 16.2 x 11.8 cm
- Nett weight: 0.312 kg

- Number of consumer packagings: 24
- Tare weight: 0.09 kg

Outer Carton

- Gross weight: 2.038 kg
- GTIN: 1 48 95229 10055 5
- Outer carton (L x W x H): 34.9 x 32.5 x 26.4 cm
- Nett weight: 1.248 kg
- Number of consumer packagings: 96
- Tare weight: 0.79 kg

Packaging dimensions

- Packaging dimensions (W x H x D):
5 x 17.2 x 2.5 cm
- EAN: 48 95229 10055 8
- Gross weight: 0.029 kg
- Nett weight: 0.013 kg
- Number of products included: 1
- Tare weight: 0.016 kg

Highlights

Neodymium speaker drivers

Neodymium is the best material for producing a strong magnetic field for greater sensitivity in a voice coil, better bass response and a pure balanced sound quality.

14.2mm speaker drivers



Quality 14.2mm speaker drivers with neodymium magnet delivers rich bass and clear sound.

Ergonomic design

This ergonomic design will fit perfectly into your ear canal

Designed for comfort

The design is based on ear geometry for comfortable wear and snug fit for everyone.

