

# PHILIPS

Replacement Foil

Bodygroom  
replacement foil


Fits BodyGroom S3000 Series

Fits BodyGroom S5000 Series

Fits BodyGroom S7000 Series

Fits Click&Style (S500/700)

TT2000/43



## Replacement Trim & Shave Foil

Change foil every 12 months for best results

Replacement shaving foil head TT2000/51 for Philips Bodygroom shaver range 3000, 5000 and 7000

**Best performance that feels great**

- Safe and less irritation for most comfortable body grooming

**Brings maximum ease of use**

- 100% waterproof for shower use and easy cleaning

# Highlights

## 100% waterproof



100% waterproof for easy use and cleaning.

## Safe and less irritation



The shaver head features patented rounded tips and a hypo-allergenic foil to protect your skin while shaving. The bi-directional trimmers cut longer hairs, which are shaved by the foil for a closer result.

# Specifications

## Cleaning

**Usage:** Use cleaning spray HQ110

## Content

**Packaging:** 1 shaving head foil

## Compatibility

**Fits product type:** QG33XX, TT20XX, BG70XX, BG30XX, BG50XX, SBG315, BG20XX, YS52X, YS53X

