

PHILIPS

Replacement Foil

Bodygroom
replacement foil


Fits BodyGroom S3000 Series

Fits BodyGroom S5000 Series

Fits BodyGroom S7000 Series

Fits Click&Style (S500/700)

TT2000/43



Replacement Trim & Shave Foil

Change foil every 12 months for best results

Replacement shaving foil head TT2000/51 for Philips Bodygroom shaver range 3000, 5000 and 7000

Best performance that feels great

- Safe and less irritation for most comfortable body grooming

Brings maximum ease of use

- 100% waterproof for shower use and easy cleaning

Highlights

100% waterproof



100% waterproof for easy use and cleaning.

Safe and less irritation



The shaver head features patented rounded tips and a hypo-allergenic foil to protect your skin while shaving. The bi-directional trimmers cut longer hairs, which are shaved by the foil for a closer result.

Specifications

Cleaning

Usage: Use cleaning spray HQ110

Content

Packaging: 1 shaving head foil

Compatibility

Fits product type: QG33XX, TT20XX, BG70XX, BG30XX, BG50XX, SBG315, BG20XX, YS52X, YS53X

