


PHILIPS

Replacement
shaving foil head

TT2000



Perfect trim and shave

Change foil every 12 months for best results

Replacement shaving foil head TT2000/43 for Philips Bodygroom shaver range TT2021 to TT2030.

Best performance that feels great

- Safe and less irritation for most comfortable body grooming

Brings maximum ease of use

- Wet and dry use; for use in shower and easy cleaning
- 100% waterproof for shower use and easy cleaning


Highlights

100% waterproof



100% waterproof for easy use and cleaning.

Safe and less irritation



Safe and less irritation for most comfortable body grooming

Wet and dry use

100% water proof Philips shaver makes it convenient to trim and shave your body in the shower and easy to clean.

Specifications

Cleaning

Usage: Use cleaning spray HQ110 for at least once a week for optimum results.

Content

Packaging: 1 shaving head foil

Replaceable part

Fits product types: TT2020/15, TT2020/30, TT2020/31, TT2020/35, TT2021/30, TT2021/32, TT2021/33, TT2021/34, TT2021/35, TT2021/36, TT2021/50, TT2022/30, QG3280/41, TT2030/36, QG3280/32, TT2030/60, TT2030/61, TT2040/32

